

Material List

Kind Of Material	Manufacture	Composition & Note	CE No.
NP-Alloy	Wirobond® 280 (Bego,Germany)	Co60.2 %, Cr25%,W6,2%,Mo4,8%,Ga2,9%,Si<1%,Mn<1%	#0197
	Wiron light (Bego,Germany)	Ni 65%,Cr 22,5%,Mo 9,5%,Si <1%,Nb <1%,Fe <1%,Ce <1%	#0197
	Co-Cr (Argen ,USA)	Co61.00% ,Cr27.00% ,Mo6.00% ,W5.00% ,Si1.00% ,Mn < 1,Fe < 1,C < 1	#0086
	Ni-Cr (Argen ,USA)	Ni(61.2%), Cr(25.8%),W9-11%,Mo(11%)	#0086
Precious Alloy	Argedent 74 (Argen,USA)	Au74%,Pt4%,Pd7.6%,Ag10.8%,Sn1%,In2.5%	#0086
	Argedent EURO (Argen, USA)	Au(40%),Pd(39.4%),Ag(10%)	#0086
	Argedent 52 (Argen,USA)	Au52.5%,Pd26.9%,Ag16%,Sn2%,In2.5%	#0086
	W-1 (Ivoclar, USA)	Pd53.3%,Ag37.7%,Si8.5%,In<1%	#0120
Metal Framework	Tritan®pure Titanium (Dentaurum,Germany)	min 99.5% titanium	#0129
	Co-Cr (Argen,USA)	Co64%,Cr28%,Mo6%,Si<1%,Mn<1%,C<1% ,Fe<1%	#0086
	Vitallium 2000 (Dentsply,USA)	Co63.1%,Cr28.5%,Mo6%	#0086
Ceramic Powder	Ivoclar Vivadent AG Liechtenstein	IPS e max Ceram® Ceramic	#0123
	Dentsply, USA	Ceramco®3 Ceramic	#0086
	Vita, Germany	Vita VMK 95® Ceramic	#0124
Composite	Shofu Dental, Japan	Ceramage® Light Activated Micro Ceramic	#0120
Acrylic Teeth	Yamahachi Dental, Japan	New ACE™ Acrylic Resin teeth	#0120
	Ivoclar Liechtenstein	Gnathostar® Acrylic Resin Teeth	#0123
Acrylic Resin	DENTSPLY (USA)	Lucitone 199	#0124
	TCS (USA)	TCS Unbreakable	#0470
Zirconia Blank	UPCERA (CHINA)	UPzir Blank	#0197
	WIELAND (Germany)	Zeno Blank	#0483